


Certificate

Polymer Technology Systems Inc.

using

A1CNow+

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.



Total Error =	8.4 mmol/mol
Bias =	-3.5 mmol/mol
Imprecision =	4.9 %
Grade =	Participated

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61:752-59)

Date of Certification: 01 January 2018

Date of Expiry: 01 January 2019



 IFCC Network Coordinator

 Dr. C.W. Weykamp