

Customer Bulletin – January 2017

Revised Range Cards: PTS Panels® Multi-Chemistry Controls (Lot Numbers MC23, MC24, MC25, and MC26) (REF 721)

Summary:

New range cards has been issued for PTS Panels Multi-Chemistry Controls. (LOT Numbers MC23, MC24, MC25, and MC26). No other lots are affected.

Single Analyte Glucose Range Update (including PTS Panels® eGLU test strips):

The new range card updates the range for glucose when utilizing single-analyte PTS Panels test strips. The updated ranges are provided below, with the change highlighted. Updated range cards are also available at www.ptsdiagnostics.com. The previous, obsoleted range cards are also provided for reference.

PTS Panels Multi-Chemistry Controls – Lot Number MC23

PTS Panels® Multi-Chemistry Control* Ranges	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #: MC23 Expiration Date: 2017 06 30				
Single-Analyte Test Strips				
Cholesterol	120-220	3.11-5.70	130-300	3.37-7.77
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU™	40-120	2.22-6.66	170-370	9.44-20.54
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	150-310	3.89-8.03
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	120-340	6.66-18.87
*HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). For assistance contact PTS Customer Service @ 877-870-5610 or Fax @ 317-870-5608. Outside of US @ 317-870-5610. Email: inforequest@cardiochek.com				

PS-002519 Rev. 9

PTS Panels Multi-Chemistry Controls – Lot Number MC24

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC24 Expiration Date: 2017 07 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-300	3.37-7.77
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	35-85	0.91-2.20
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	40-120	2.22-6.66	170-370	9.44-20.54
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-50	0.39-1.30	30-85	0.78-2.20
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	120-340	6.66-18.87
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). ** For CardioChek Plus analyzers version 1.09 and above. PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 10

PTS Panels Multi-Chemistry Controls – Lot Number MC25

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC25 Expiration Date: 2018 03 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-315	3.37-8.16
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	30-80	0.78-2.07
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	40-120	2.22-6.66	170-370	9.44-20.54
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	115-225	2.98-5.83	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-55	0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	110-330	6.11-18.32
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). ** For CardioChek Plus analyzers version 1.09 and above. PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 10

PTS Panels Multi-Chemistry Controls – Lot Number MC26

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC26 Expiration Date: 2018 07 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-315	3.37-8.16
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	30-80	0.78-2.07
Triglycerides	100-195	1.13-2.20	160-270	1.81-3.05
Glucose/eGLU	40-120	2.22-6.66	170-370	9.44-20.54
Multi-Analyte Test Strips				
Cholesterol	115-225	2.98-5.83	150-310	3.89-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-55	0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	150-330	1.70-3.73
Glucose	20-90	1.11-5.00	110-330	6.11-18.32
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 11

Previous Range Card Information (Obsolete):

PTS Panels Multi-Chemistry Controls - Lot Number MC23 - OBSOLETE

PTS Panels® Multi-Chemistry Control* Ranges	Level 1		Level 2	
	mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #: MC23 Expiration Date: 2017 06 30				
Single-Analyte Test Strips				
Cholesterol	120-220	3.11-5.70	130-300	3.37-7.77
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU™	50-110	2.78-6.11	190-320	10.55-17.76
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	150-310	3.89-8.03
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	120-340	6.66-18.87
<p>*HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). For assistance contact PTS Customer Service @ 877-870-5610 or Fax @ 317-870-5608. Outside of US @ 317-870-5610. Email: inforequest@cardiochek.com</p>				
				PS-002519 Rev. 9

PTS Panels Multi-Chemistry Controls - Lot Number MC24 - OBSOLETE

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC24 Expiration Date: 2017 07 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-300	3.37-7.77
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	35-85	0.91-2.20
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	50-110	2.78-6.11	190-320	10.55-17.76
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-50	0.39-1.30	30-85	0.78-2.20
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	120-340	6.66-18.87
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). ** For CardioChek Plus analyzers version 1.09 and above. PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 10

PTS Panels Multi-Chemistry Controls - Lot Number MC25 - OBSOLETE

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC25 Expiration Date: 2018 03 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-315	3.37-8.16
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	30-80	0.78-2.07
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	45-110	2.50-6.11	190-320	10.55-17.76
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77
Multi-Analyte Test Strips				
Cholesterol	115-225	2.98-5.83	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-55	0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	110-330	6.11-18.32
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). ** For CardioChek Plus analyzers version 1.09 and above.				
PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 10

PTS Panels Multi-Chemistry Controls - Lot Number MC26 - OBSOLETE

pts panels[®] Multi-Chemistry Controls Range Card	Level 1		Level 2	
	mg/dL	mmol/L	mg/dL	mmol/L
Lot #: MC26 Expiration Date: 2018 07 31				
Single-Analyte Test Strips				
Cholesterol	110-220	2.85-5.70	130-315	3.37-8.16
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	30-80	0.78-2.07
Triglycerides	100-195	1.13-2.20	160-270	1.81-3.05
Glucose/eGLU	45-110	2.50-6.11	190-330	10.55-18.32
Multi-Analyte Test Strips				
Cholesterol	115-225	2.98-5.83	150-310	3.89-8.03
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-55	0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	150-330	1.70-3.73
Glucose	20-90	1.11-5.00	110-330	6.11-18.32
* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722). PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com				

PS-002519 Rev. 11

Questions/Comments?

Please call Customer Service at +1-317-870-5610 (direct) or 1-877-870-5610 (toll-free US only).

CB 17-004 r0 01/17