

pts panels[®]

HDL Cholesterol Controls Range Card

		Level 1		Level 2	
Lot #:	HC25	mg/dL	mmol/L	mg/dL	mmol/L
Expiration Date: 2019 06 30					
Single Analyte Test Strips					
HDL Cholesterol* <i>(CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below)</i>		≤25-45	≤0.65-1.17	30-60	0.78-1.55
Multi-Analyte Test Strips					
HDL Cholesterol* <i>(CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below)</i>		≤15-50	≤0.39-1.30	30-70	0.78-1.81

* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610
Email: customerservice@ptsdiagnostics.com

PS-002516 Rev. 7