

pts panels[®]**HDL Cholesterol Controls Range Card**Lot #: **HC26**Expiration Date: **2019 08 31****Level 1****Level 2**

mg/dL

mmol/L

mg/dL

mmol/L

Single Analyte Test StripsHDL Cholesterol*
(*CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below*)

≤25-45

≤0.65-1.17

30-70

0.78-1.81

Multi-Analyte Test StripsHDL Cholesterol*
(*CardioChek[®] PA analyzers and CardioChek[®] Plus analyzers version 1.08 or below*)

≤15-50

≤0.39-1.30

30-70

0.78-1.81

* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610
Email: customerservice@ptsdiagnostics.com

PS-002516 Rev. 7