

# pts panels<sup>®</sup>

HDL Cholesterol Controls **Range Card**

**Lot #:** HC27

**Expiration Date:** 2020-04-30

## Single Analyte Test Strips

HDL Cholesterol\*  
(*CardioChek<sup>®</sup> PA analyzers and CardioChek<sup>®</sup> Plus analyzers version 1.08 or below*)

## Multi-Analyte Test Strips

HDL Cholesterol\*  
(*CardioChek<sup>®</sup> PA analyzers and CardioChek<sup>®</sup> Plus analyzers version 1.08 or below*)

## Level 1

mg/dL

mmol/L

≤25-45

≤0.65-1.17

≤15-50

≤0.39-1.30

## Level 2

mg/dL

mmol/L

30-70

0.78-1.81

30-70

0.78-1.81

\* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610  
Email: customerservice@ptsdiagnostics.com

PS-002516 Rev. 7