

**pts panels<sup>®</sup>**

Multi-Chemistry Controls      Range Card

Level 1

Level 2

Lot #:                    MC24

Expiration Date:        2017 07 31

mg/dL

mmoL/L

mg/dL

mmoL/L

**Single-Analyte Test Strips**

Cholesterol	110-220	2.85-5.70	130-300	3.37-7.77
HDL Cholesterol** <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	35-85	0.91-2.20
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	40-120	2.22-6.66	170-370	9.44-20.54
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77

**Multi-Analyte Test Strips**

Cholesterol	110-220	2.85-5.70	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>	15-50	0.39-1.30	30-85	0.78-2.20
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	120-340	6.66-18.87

\* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).

\*\* For CardioChek Plus analyzers version 1.09 and above.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 10