

pts panels[®]

Multi-Chemistry Controls Range Card

Level 1

Level 2

Lot #: MC25

Expiration Date: 2018 03 31

mg/dL

mmol/L

mg/dL

mmol/L

Single-Analyte Test Strips

Cholesterol	110-220	2.85-5.70	130-315	3.37-8.16
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-45	≤0.65-1.17	30-80	0.78-2.07
Triglycerides	90-180	1.02-2.03	150-260	1.70-2.94
Glucose/eGLU	40-120	2.22-6.66	170-370	9.44-20.54
Ketone	≤2-20	≤0.19-1.92	20-60	1.92-5.77

Multi-Analyte Test Strips

Cholesterol	115-225	2.98-5.83	150-310	3.89-8.03
HDL Cholesterol** <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	15-55	0.39-1.42	30-90	0.78-2.33
Triglycerides	70-190	0.79-2.15	130-300	1.47-3.39
Glucose	20-90	1.11-5.00	110-330	6.11-18.32

* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).

** For CardioChek Plus analyzers version 1.09 and above.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 10