

pts panels[®] Multi-Chemistry Controls Range Card		Level 1		Level 2	
		mg/dL	mmoL/L	mg/dL	mmoL/L
Lot #:	MC31				
Expiration Date:	2020 03 31				
Single-Analyte Test Strips					
Cholesterol		120-230	3.11-5.96	150-330	3.89-8.55
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤25-50	≤0.65-1.30	30-80	0.78-2.07
Triglycerides		100-200	1.13-2.26	180-300	2.03-3.39
Glucose/eGLU		40-130	2.22-7.22	160-370	8.88-20.54
Multi-Analyte Test Strips					
Cholesterol		110-230	2.85-5.96	140-330	3.63-8.55
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤15-60	≤0.39-1.55	30-90	0.78-2.33
Triglycerides		90-200	1.02-2.26	150-320	1.70-3.62
Glucose		20-100	1.11-5.55	125-350	6.94-19.43
<p>* For CardioChek, CardioChek PA, and CardioChek Plus (version 1.08 or below) analyzers, HDL Cholesterol requires a separate HDL Cholesterol Control (Cat. No. 722).</p> <p>PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com</p>					
PS-002519 Rev. 11					